

Supplementary Table 1

Primers used for *TP53* mutation screening.

Exon	Sequence
<i>TP53</i> exon 4.1 forward	cctggtcctctgactgctct
<i>TP53</i> exon 4.1 reverse	ggtgttaggagctgctggtg
<i>TP53</i> exon 4.2 forward	gtccagatgaagctcccaga
<i>TP53</i> exon 4.2 reverse	gccaggcattgaagtctcat
<i>TP53</i> exon 5 forward	tctgtctccctcccttcata
<i>TP53</i> exon 5 reverse	aaccagccctgtcgctct
<i>TP53</i> exon 6 forward	caggcctctgattcctcact
<i>TP53</i> exon 6 reverse	cttaaccctcccccagag
<i>TP53</i> exon 7 forward	cttgggcctgtgttatctcc
<i>TP53</i> exon 7 reverse	gggtcagaggcaagcaga
<i>TP53</i> exon 8 forward	gcttcttgcctgtttgttt
<i>TP53</i> exon 8 reverse	gggagtagatggagcctgg
<i>TP53</i> exon 9 forward	ggagaccaagggtgcagttt
<i>TP53</i> exon 9 reverse	ccccaaattgcaggtaaaca